

Physician Chair Scales

Frequent weight checks of frail, handicapped, and elderly patients can be a tremendous, time-consuming task without the proper equipment. Now Detecto simplifies that job with a selection of innovatively designed chair scales. Developed to ensure patient safety and comfort as well as to reduce weighing time, Detecto chair scales are conveniently

portable and offer easy side or front access for your special patients. Choose from mechanical or digital models, each constructed as a singular, self-contained unit. Both types feature movable arms and footrest and oversized wheels for easy maneuverability.



**Mechanical
Model 475**



**Digital
Model 6475**

Arms raise easily on both models so patient can be conveniently assisted into and out of the chair scale.



Waist-high LCD readout on the digital model allows rapid weighing with accuracy to .2 lb or .1 kg.

Choose Detecto Chair Scales...reliably built for patients with special health care problems

Sturdy construction combines with time-saving design features in the most functional chair scales available anywhere. Complete with heavy-duty understructure, these scales come fully assembled and ready to use. Pleasing tones blend with any clinical environment, and these units can be used in relatively small spaces - on carpet or tile.



Oversized wheels facilitate movement of scale from place to place...especially over thresholds and elevator entries.



Lift-away arms allow patient to be positioned in chair scale with ease.



Footrest lifts up and out of the way for safe patient seating as well as for convenient transporting of scale.

Mechanical

Mechanical chair scales are available in pounds or pounds/kilos models. Die-cast, easy-to-read beams are readable from front and back. Black vinyl upholstery is complimented by the metal supporting structure.



Specifications

Model	Capacity
475	400 lb x 4 oz
4751	400 lb x 4 oz / 180 kg x 100 g

Beam: Die-cast, easy-to-read from both sides

Features: Lift-away arms and footrest
Portability
Oversized wheels for easy transport

All Detecto scales are engineered for dependable accuracy and are virtually maintenance-free.

Make Detecto your scale of choice when it comes to caring for your special patients.

All Detecto scales are tested with weights traceable to the US Bureau of Standards.

We reserve the right to improve, enhance or modify features and specifications without prior notice.

Digital

Easy-to-read digital rapidly displays weight to .2 lb or .1 kg. Battery power enhances this electronic scale's portability. 115 VAC power adapter available as an option.



Specifications

Model	Capacity
6475	400 lb x .2 oz
6475KG	180 kg x .1 kg

Display: .7" high-contrast LCD

Power: 6 "C" size Alkaline batteries (not included)

Option: 115 VAC power adapter



P.O. Box 151, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
Visit our Web site at <http://www.detectoscale.com>

Sold By: