




PRODUCT SPECIFICATIONS

INFANT SCALE MODEL 6730



Detecto has a Weigh with Babies

Designed to be used with utmost ease, Model 6730 Infant Scale can weigh even the most active baby. With step-by-step instructions right on the keypad, there's never any question about what to do next.

<p>Dimensions: 23" W x 14" D x 7" H</p> <p>Tray: 25 3/4" W x 12" D at top, tapering to 25 3/4" W x 9 1/4" D at bottom</p> <p>White molded plastic</p> <p>Capacity: 30 lb x .2 oz/15 kg x 5 g</p> <p>Power: 8 "D" size Alkaline batteries or 8 "D" size Ni-Cad batteries (batteries not included)</p> <p>Controls: On/off, lb/kg, zero (used to return to zero or as tare to offset weight of blanket or cover.)</p>	<p>Operating Temperature: 0½ to 40½ C (32½ to 104½ F)</p> <p>Display: 5 digits, 1" high LCD</p> <p>Status Indicators: Pounds/ounces, kilograms/grams, weight/lock, low battery</p> <p>Functions: Auto power off • Lock/release • Automatic zero</p> <p>Ship Weight: 25 lb</p> <p>Weight Sensor: Strain-gauge load cell</p> <p>Motion Detector: Senses motion on scale and locks in weight display after motion has ceased or after 5 seconds</p> <p>Color: White tray, white front panel with blue accenting</p>
--	--



Cardinal Scale reserves the right to improve, modify, or enhance features and specifications without prior notice.

