


BMI Bar Graph

Hold/Release patient weight



Up to 99 tares (optional)

9-digit patient ID


Patient ID Barcode Scanning

Send patient ID to 758C indicator for EMR/EHR records



EMR/EHR Connectivity

Efficiently send patient data to a PC for electronic recordkeeping via the 758C indicator's serial, USB, or optional Ethernet port or to a printer.



Digital height rod input for optimum measurement efficiency

BATTERY ACCESS
USE 6 "C" CELL BATTERIES

6 C batteries (not included) or optional AC adapter

Long-lasting, durable steel construction

