



Waist-High Digital Physician Scale

Designed and constructed to provide years of reliable service, our waist-high, digital physician scale features the at-a-glance convenience and accuracy of an LCD weight display.

Model 8437

Capacity:	400 lb x .2 lb / 180 kg x .1 kg
Platform size:	10.5" x 14.5" / 267 mm x 368 mm
Display:	.7" / 17.8 mm high-contrast LCD
Color:	White
Power:	6 "C" size alkaline batteries
Keys:	On/Off, Zero, Lb/Kg, Hold/Release Body mass index calculator

Optional height rod available



Cardinal Scale reserves the right to improve, modify, or enhance features and specifications without prior notice.

