

Model 8440 Infant Scale



*Incredibly Low Price
Extremely High Accuracy
Exceptionally Small Footprint*

The Detecto 8440 digital infant scale provides accurate, safe weighing of infants and features selectable weight units, lb-oz or kilograms, display hold, and a memory function. Compact, battery powered and portable, the 8440 is perfect for visiting nurses, home or wherever space is limited and performance is paramount.

Large, easy to read display with lock-on weight!

Ergonomically designed for both ease of use and infant safety.



All the displays that you need!



Displays weight in lb/oz or kg at the touch of a button.



Built-in clock shows time and date for recording purposes.



Detachable tray for easy cleaning and portability.



Specifications and Features

- | | |
|-------------------|--|
| Capacity: | 44 lb x 1/2 oz / 20 kg x 10 g |
| Display: | 4-3/4 digits LCD 1.0" (25mm) high |
| Clock function: | Time and date |
| Memory function: | Memory for comparison of current and stored weight |
| Accuracy: | 0, 2% +/- 1 digit of reading over 200 grams |
| Power Source: | 6V DC (4 - AA batteries or included AC adapter) |
| Tray dimensions: | 20 1/2" L x 10 3/4" W
(520mm L x 272mm W) |
| Scale dimensions: | 13" L x 9 3/4" W x 7 1/4" H
(333mm L x 246mm W x 184mm H) |
| Weight: | 8 lb (3.5 kg) |

Cardinal Scale reserves the right to improve, enhance or modify features and specifications without prior notice.

 **DETECTO**[®]
A Division of Cardinal Scale Manufacturing Co.

P.O. Box 151, Webb City, MO 64870 USA

Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001

Visit our Web site at www.detectoscale.com

© Copyright 2002 Cardinal Scale Mfg. Co. • Printed in USA • BLC/5M/0502/C-164

Sold By: