



model **8440**
BABY
SCALES



*Detecto
Model 8440
Baby Scale*

Providing accurate and safe weighing of infants and toddlers, Detecto's model 8440 digital scale is ideal for home health and medical clinics. The 8440's weighing tray may be removed for standing toddler weighing and for easier transport. The 1-in / 25-mm high LCD offers pounds and kilograms readouts plus features a built-in clock. The scale's Weight Difference key provides the weight difference between the current weight on the scale and a saved weight in memory along with its time and date.

*Detecto
Model
8440-CASE
Carrying
Case*



CAPACITY: 44 lb X 1/2 oz / 20kg x 0.01 kg

8.9 1/2 lb. oz

DETECTO
www.Detecto.com
MODEL 8440 BABY SCALE

1-in / 25-mm
high bold LCD
with fractions

Ergonomically
designed for both
ease of use and
infant safety

TARE

Large, easy-to-read display with weight recall



The tray collapses for easier transport and to allow for convenient standing toddler weighing.

- Accuracy to 1/2 ounce
- Storage of previous weight memory recall
- Calculation of change in weight
- Low battery indicator
- Optional carrying case available (8440-CASE)

8440 BABY SCALE SPECIFICATIONS

CAPACITY: 44 lb x 1/2 oz / 20 kg x 10 g

DISPLAY: 1 in / 25 mm high LCD

TRAY DIMENSIONS: 21 in W x 9.5 in D x 3.5 in H /
53.3 cm W x 24.1 cm D x 8.9 cm H

SCALE DIMENSIONS: 13.7 in W x 11 in D x 6.9 in H /
35 cm W x 28 cm D x 18 cm H

SHIPPING WEIGHT: 8 lb / 3.6 kg

POWER: (1) 9V Alkaline battery (not included) and (1) 3V
CR2032 lithium back-up battery (included) for memory

CLOCK FUNCTION: Time and date

MEMORY FUNCTION: Memory for comparison of current and stored weight



Selectable weight units (lb or kg)



Date and time displays automatically
in sleep mode or by pressing the
Time/Weight button during weighing

Detecto reserves the right to improve, enhance, or modify features and specifications without prior notice.

 **DETECTO**
A Division of Cardinal Scale Manufacturing Co.

203 E. Daugherty, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
www.Detecto.com

Sold By: