

**PUT YOUR CONFIDENCE IN...  
AND BABIES ON...**

# **Detecto's Model 8450 Baby Scale**



**Large, Easy-to-read, 1" LCD Display  
Built-In Measuring Tape  
Battery or AC Power Adapter  
Memory and Lock-in Weight Features  
Light Weight, One Piece Construction  
Accurate Weights Displayed in lb/oz or kg**

# DETECTO'S 8450 BABY SCALE...

## Perfect for Hospitals, Physicians' Offices or, at Home

This extremely light weight scale can be used either in a permanent location or as a portable scale. The platform, base and padded carrying handles are incorporated into one complete unit making it easy to transport from place to place. Padded side rails and handle covers are slightly contoured ensuring baby will be cradled comfortably & securely while being weighed. Measure infant's length in inches or centimeters with built-in measuring tape.



The 8450 is extremely easy-to-use. Accurate weight is displayed in pounds/ounces or kilograms on the large, 1" LCD. Conversion is as simple as the touch of a button. Optional units, pounds only and ounces, are available via setup. The Memory feature allows operator to record weight. Lock-In weight display means baby's movements won't effect accurate weight readings.

| Specifications         |   |
|------------------------|---|
| <b>Capacities:</b>     | 20 kg x 0.01 kg/44 lb x 0.02 lb/704 oz x 0.5 oz |
| <b>Weight Units:</b>   | lb/oz, kg, lb, oz                               |
| <b>Power:</b>          | 9V battery or 9V/100mA AC adapter               |
| <b>Display:</b>        | 1" LCD display                                  |
| <b>Operation</b>       |   |
| <b>Temperature:</b>    | 5° to 35°C/41° to 95°F                          |
| <b>Keys:</b>           | On/Zero/Off, kg, lb, oz, lb/oz, hold            |
| <b>Measuring Tape:</b> | Inch and centimeter graduations                 |
| <b>Platform:</b>       | 500 x 350 x 50 mm/20" x 14" x 2"                |
| <b>Overall</b>         |   |
| <b>Dimensions:</b>     | 620 x 350 x 50 mm/24" x 14" x 2"                |
| <b>Net Weight:</b>     | 3.5 kg  |
| <b>8450-CASE:</b>      | optional carrying case                          |

We reserve the right to improve, enhance or modify features and specifications without prior notice.

 **DETECTO**  
 A Division of Cardinal Scale Manufacturing Co.  
 P.O. Box 151, Webb City, MO 64870 USA  
 Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001  
 Visit our Web site at [www.detectoscale.com](http://www.detectoscale.com)

**Sold By:**