

DR550C

Stainless Steel Portable Floor Scale



DR550C SPECIFICATIONS

Capacity:	550 lb x 0.2 lb / 250 kg x 0.1 kg
Platform Size:	12.2" W x 11.8" D x 1.97" H / 31 cm W x 30 cm D x 5 cm H
Display:	4 digit, 7 segment, liquid crystal display, 1"/25mm high
Keypad:	4 keys - LB/KG, TARE/HOLD, ON, OFF
Platform Construction:	Stainless steel
Power:	6 "AA" batteries (not included) or 115 volt AC adapter (included)
Shipping Weight:	13 lb / 6 kg

Accurate and reliable, Detecto's economical DR550C stainless steel platform scale is light-weight and portable, making it perfect for mobile clinics and home care nurses. The remote indicator features a large 1"/25mm high LCD display, units conversion, and tare. The unit has a slip-resistant mat to ensure patient safety when getting on and off the scale. The DR550C is battery powered, so you can take it anywhere you need it.

Optional Black Carrying Case

Model DR400C-CASE

Carrying case with handle protects the DR550C scale for easy transport.

