

HOME HEALTH CARE SCALES



EYE-LEVEL PHYSICIAN SCALE

Model 437

Capacity: 400 lb x 4 oz
Platform Size: 10.5" W x 14.5" D / 27 cm W x 37 cm D
Scale Height: 59" / 150 cm



MECHANICAL DIAL PERSONAL FITNESS SCALE

Model 080

Capacity: 300 lb x 1 lb
Platform Size: 11" W x 13" D / 28 cm W x 33 cm D
Scale Height: 36" / 91 cm

DIGITAL BABY AND TODDLER SCALE

Model 8440

Capacity: 44 lb x .5 oz / 20 kg x 10 g
Tray Size: 20.5" L x 10.75" W / 52 cm L x 27 cm W
Platform Size: 9.5" L x 6.5" W / 24 cm L x 17 cm W



VISITING NURSE SCALE

Model DR400C

Capacity: 400 lb x 0.5 lb / 181 kg x 0.2 kg
Platform Size: 12" W x 12" D / 31 cm W x 31 cm D



PERSONAL FLOOR SCALE

Model D315

Capacity: 330 lb x 1 lb / 150 kg x 500 g
Dimensions: 13" W x 18.5" L x 4.75" H /
33 cm W x 47 cm L x 12 cm H



PERSONAL FLOOR SCALE

Model D1130

Capacity: 300 lb x 1 lb
Dimensions: 11.5" W x 17.25" L x 3" H /
29 cm W x 44 cm L x 8 cm H



DETECTO[®]
A Division of Cardinal Scale Manufacturing Co.