

# LONG-TERM CARE SCALES



## FOLDING PORTABLE WHEELCHAIR SCALE

**Model 6550**

Capacity: 800 lb x 0.2 lb / 360 kg x 0.1 kg  
Platform Size: 28" W x 32" D / 71 cm W x 81 cm D



## ROLLING CHAIR SCALE

**Model 6880**

Capacity: 440 lb x 0.2 lb / 200 kg x 0.1 kg  
Seat Size: 17.5" H x 17.5" W / 44 cm H x 44 cm W  
Zero Turning Radius Wheels



## GERIATRIC WHEELCHAIR SCALE

**Model FHD-133-II**


Capacity: 1,000 lb x 0.2 lb / 450 kg x 0.1 kg  
Platform Size: 36" W x 36" D /  
91 cm W x 91 cm D



## "ROLL-A-WEIGH" WHEELCHAIR SCALE

**Model CR1000D**

Capacity: 1,000 lb x 0.5 lb / 450 kg x 0.2 kg  
Platform Size: 30" W x 32.5" D /  
76 cm W x 83 cm D



## FLIP SEAT CHAIR SCALE

**Model 6876**

Capacity: 600 lb x 0.2 lb / 270 kg x 0.1 kg  
Platform Size: 24" W x 24" D / 61 cm W x 61 cm D



## HEAVY-CAPACITY IN-BED SCALE

**Model IB600**

Capacity: 600 lb x 0.2 lb / 272 kg x 0.1 kg  
Stretcher: 6' / 1.8 m long, Thermoplastic Sure-Chek® fabric



## DIGITAL CHAIR SCALE

**Model 6475**

Capacity: 400 lb x 0.2 lb / 180 kg x 0.1 kg  
Seat Size: 17.5" H x 17.5" W / 44 cm H x 44 cm W



**DETECTO®**  
A Division of Cardinal Scale Manufacturing Co.