

SPECIALTY SCALES FOR SURGERY



HIGH-CAPACITY STAND-ON SCALES

Model 6856

Capacity: 800 lb x 0.2 lb /
360 kg x 0.1 kg
Platform Size: 24" W x 24" D /
61 cm W x 61 cm D

Model 6855

Capacity: 600 lb x 0.2 lb /
270 kg x 0.1 kg
Platform Size: 18" W x 14" D /
46 cm W x 26 cm D



BARIATRIC FLIP SEAT SCALE

Model 6868

Capacity: 800 lb x 0.2 lb / 360 kg x 0.1 kg
Platform Size: 32" W x 34" D / 81 cm W x 86 cm D



DIAPER/SPONGE SCALE

Model AP-4KD

Capacity: 3,999 g x 1 g
Display: 4 digit, 0.5"/13 mm high red LED
Overall Size: 2.75" H x 7.5" W x 7" D /
7 cm H x 19 cm W x 18 cm D
Tray Size: 2" H x 9" W x 6" D /
5 cm H x 23 cm W x 15 cm D

Construction: Stainless steel scale with stainless steel tray



DIAPER/SPONGE SCALE

Model PSS-2KD

Capacity: 70 oz x 0.1 oz / 1,998 g x 2 g
Display: 4 digit, 0.5"/13 mm high LCD
Overall Size: 2.75" H x 7.5" W x 7" D /
7 cm H x 19 cm W x 18 cm D
Tray Size: 2" H x 9" W x 6" D /
5 cm H x 23 cm W x 15 cm D

Construction: Scale is metal with baked enamel finish and stainless steel tray



TOP LOADING DIAL SCALE

Model PT-5R

Capacity: 5 lb x 1/2 oz
Dial Size: 6" / 15 cm diameter
Overall Size: 7.5" H x 5.25" W x 7.25" D /
19 cm H x 13 cm W x 18 cm D
Platform: 7.875" W x 5.875" D /
20 cm W x 15 cm D stainless steel



A Division of Cardinal Scale Manufacturing Co.