

# VETERINARY PRODUCTS


## HIGH CAPACITY DIGITAL VET SCALE

### Model VET400

Capacity: 400 lb x 0.2 lb / 200 kg x 0.1 kg

Platform Size: 1.75" H x 42" W x 22" D

Power: 6 "C" size Alkaline or Nickel-Cadmium batteries or optional AC adapter



## MEDIUM CAPACITY DIGITAL VET SCALE

### Model VET50

Capacity: 50 lb x 0.02 lb / 20 kg x 0.01 kg

Platform Size: 3" H x 14" W x 12" D

Power: 6 "C" size Alkaline or Nickel-Cadmium batteries or optional AC adapter



## BAKED EPOXY STEEL STEP-ON CANS

### Model P

Capacities: From 16 quart (4 gallon)  
to 100 quart (25 gallon)

Dimensions: From 13" H x 11.75" W x 13" D  
to 27.5" H x 16.75" W x 17.75" D



## UNIVERSAL GLOVE BOX HOLDERS

### Model GH

Capacity: Models available that hold 1, 2, and 3 boxes

Dimensions: From 10.5" H x 6.5" W x 4.25" D  
to 10" H x 18" W x 4" D



A Division of Cardinal Scale Manufacturing Co.