

475 / 4751

Weigh Beam Chair Scale

Frequent weight checks of frail, handicapped, and elderly patients can be a tremendous, time-consuming task without the proper equipment. Detecto simplifies that job with a selection of innovatively-designed chair scales. Developed to ensure patient safety and

comfort as well as to reduce weighing time, Detecto chair scales are conveniently portable and offer easy side or front access for your special patients. The arms raise easily on both models, so the patient can be conveniently assisted into and out of the chair scale.



*Die-cast, easy-to-read
beams are readable
from the front and back.*

SPECIFICATIONS

Capacity:

475 400 lb x 4 oz
4751 400 lb x 4 oz /
175 kg x 100 g

Seat Size:

17.5" H x 17.5" W /
44 cm H x 44 cm W

Shipping Weight:

78 lb / 35 kg