

Detecto Medical Scale Indicators



MODEL 750



MODEL 758C

DETECTO manufactures the most comprehensive line of USA-made medical scale products in the clinical industry. The models 758C and 750 weight displays provide you with versatility to meet the demands of each unique weighing application.

The model 750 display features an easy-to-use keypad and extra-large display perfect for patient self-weighing and viewing by the attending nurse or physician. A quick and simple BMI calculator with a scroll feature to input height allows you the ease of BMI calculation when you require it. RS232 serial port and easy-to-replace "AA" batteries makes the 750 ideal for clinics, hospitals, and physician offices.

From wheelchair-based applications to electronic medical record updates, the 758C is the most versatile clinical scale weight display available. A full numeric keypad allows fast entry of wheelchair weight. For those applications requiring weight, height and BMI percentages simultaneously, the 758C is the right choice. The 758C also includes EMR/EHR-ready RS232 serial, USB or optional Ethernet ports.



750 Indicator

- Extra-large, 0.75 inch high easy-to-read LCD display.
- RS232 serial output for EMR/EHR connectivity.
- Just touch the BMI key and utilize the toggle arrows to adjust the patient height for the BMI calculation.
- Compact design is ergonomic and easy for the patient to weigh.
- Use either 6 "AA" batteries or the AC adapter (optional).
- Easy to remove from the scale column for battery replacement.
- Lock/Release key for viewing weight after the patient steps away from the scale.
- Units conversion for LB or KG.
- Steel construction for years of long-lasting service.
- Made in USA for convenient parts and service requirements.



758C Indicator

- 0.5 Inch high graphic display for viewing of weight, height and BMI percentage calculations simultaneously.
- USB, RS232 serial, and optional Ethernet ports for EMR/EHR connectivity. Custom programming available to connect with your EMR software package (consult factory).
- 11-digit patient ID entry for BMI tracking.
- Full numeric keypad for entry of tare weight (weight of empty wheelchair).
- Lock/Hold feature for viewing weight after the patient steps from the scale.
- Choice of 6 "C" sized batteries and AC adapter (optional).
- Steel construction for years of long-lasting service.
- Made in USA for parts replacement and factory support.

