

EU CLASS III DIGITAL SCALES



PHYSICIAN SCALE MODEL EU522

Capacity: 200 kg x 0.1 kg
Platform: 26 cm W x 34 cm L
Scale Height: 95 cm
*Lock/Release key for viewing weight after the patient steps away from the scale.
Body Mass Index calculation.*

CHAIR SCALE MODEL 6875KGEU

Capacity: 200 kg x 0.1 kg
*Lift-up arms and footrest.
Contoured, sloping back rest allows the patient to assume a relaxed position while weighing.*



SCALE WITH HEIGHT ROD MODEL EU522HR

Capacity: 200 kg x 0.1 kg
Platform: 26 cm W x 34 cm L
Scale Height: 95 cm
Height Rod: Measures 76 cm - 200 cm
Slip-resistant base cover for patient safety and assurance.



PORTABLE SCALE MODEL 6876KGEU

Capacity: 270 kg x 0.1 kg
Platform: 61 cm x 61 cm
The seat and armrests both flip up to allow weighing sitting or standing. Built-in wheels for easy transport.



FLOOR SCALE MODEL EU544

Capacity: 200 kg x 0.1 kg
Platform: 26 cm W x 34 cm L
*Sleep Mode - Auto ON/OFF to preserve battery life.
Slip-resistant base cover.*



PORTABLE SCALE MODEL 6880KGEU

Capacity: 200 kg x 0.1 kg
*Lift-up arms and footrest.
Zero turning radius ability with built-in wheels for tight turns.*





EU CLASS III DIGITAL SCALES



FOLD-UP WHEELCHAIR SCALE MODEL 6550KGEU

Capacity: 300 kg x 0.1 kg

Platform: 76 cm W x 81 cm D (flat portion)

The vertical fold-up design requires minimal floor space and allows easy transportation with dual direction wheels.



BARIATRIC SCALE MODEL 6868KGEU

Capacity: 300 kg x 0.1 kg

Platform: 81 cm W x
86 cm L

Unique padded flip-up seat allows weighing in a sitting or standing position.



WHEELCHAIR SCALE MODEL 6500KGEU

Capacity: 300 kg x 0.1 kg

Platform: 76 cm W
x 81 cm D (flat portion)

Built-in ramps for entering and exiting the scale. Integral wheels and carrying handles.



BARIATRIC SCALE MODEL 6855KGEU

Capacity: 300 kg x 0.1 kg

Platform: 46 cm W x 26 cm L

This waist-high digital medical scale features built-in wheels for easy transport and padded handgrips for patient stability.



PORTABLE SCALE MODEL CR-1000KGEU

Capacity: 450 kg x 0.2 kg

Platform: 76 cm W x
83 cm L

Features built-in wheels, carrying handles, wheelstop, and lightweight design for easy transport.



BARIATRIC SCALE MODEL 6856KGEU

Capacity: 300 kg x 0.1 kg

Platform: 61 cm W x 61 cm L

The sturdy stainless steel handrails and padded foam handgrips provide patient assurance and stability in weighing.




WHEELCHAIR SCALE MODEL BRW1000KGEU

Capacity: 450 kg x 0.2 kg

Platform: 102 cm W x 76 cm L

This portable EMR-ready digital bariatric scale accommodates new extra-wide wheelchairs.



DETECTO[®]
A Division of Cardinal Scale Manufacturing Co.

DETECTO EU WAREHOUSE

Unit 17/18 Mercers Road, Chapel Pond Hill,
Bury St. Edmunds, Suffolk, IP32 7HX, England

TEL: +44 1284 703117 • FAX: +44 1284 703559

E-MAIL: sales@detectoscale.co.uk • WEB: www.Detecto.com