

GP-400-758C

Portable Low-Profile Healthcare Scale



Ideal for home healthcare, uses for these lightweight, portable scales are almost limitless. With detachable indicators and built-in carrying handles, they go with you anywhere weighing is needed. Their large, easy-to-clean platforms (2" H x 18" W x 14" D / 6 cm H x 46 cm W x 36 cm D) feature a black slip-proof mat and leveling eye. **MADE IN USA. LEGAL FOR TRADE.**

MODEL	CAPACITY	SHIP WT.
GP-400-758C	400 lb x 0.2 lb / 200 kg x 0.1 kg	31 lb / 14 kg
6800-1045	AC Adapter	1 lb / 453 g

EMR-READY

The scale's made-in-the-USA weight indicator includes DB9 serial port, USB, or optional Ethernet output for data transfer to a printer or a PC-based EMR (electronic medical record) software program. With the time/date feature, the scale can send the weight, time, date, BMI (Body Mass Index), and height directly to a printer or a PC.



Certificate No.
758C Display: 06-104
GP-400 Base: 94-098

