

For Wheelchair,
Chair, Stretcher,
and Stand-On
Weighing



STREAMLINE YOUR WEIGHING




DETECTO.
www.Detecto.com

DETECTO In-Floor Medical Scales



- Extra-Wide Platform to Accomodate Wheelchairs, Dialysis Chairs, Beds, and Ambulatory Stretchers.
- Pit Frame Provided.
- Flush or Wall-Mounted 1" High LCD Readouts.
- 1,000 lb x 0.2 lb / 455 kg x 0.1 kg Capacity.
- Multiple Platform Sizes Available.

DETECTO'S IN-FLOOR SCALES ARE EXCEPTIONAL FOR DIALYSIS PATIENT MONITORING AND BARIATRIC WEIGHING. The platform is level with the floor for easy access, so it also works well for multi-purpose weighing. The patient simply walks or is wheeled onto the platform and weighed. Stretchers and beds fit easily onto the 6-foot long model with no need for dual scale systems. The digital indicator is AC powered for continual use and gives you a choice of pound or kilogram weight readings. Handrail optional.



Digital readout with paper tape printer

DETECTO MANUFACTURES A FULL LINE OF MEDICAL PRODUCTS



Bariatric Scales



Wheelchair Scales



Physician Scales



Glove Box Holders



Step-On Trash Cans



In-Bed Scales

