

Pediatric & Infant Scales




**Model
8435**



**Model
8432-CH**



**Model
459-CH**



**Model
450**



**Model
243**



**Model
253**



In a hospital, pediatrician's office, clinic, or the home nursery, there is a style to fit your need. Detecto's infant and pediatric scales are available in a choice of pounds or kilograms - mechanical, 110 VAC, or battery operation. All are engineered with the active baby in mind and built with the uncompromising quality you expect from Detecto.





DETECTO'S MECHANICAL and DIGITAL PEDIATRIC and INFANT SCALES

450 Series mechanical baby scales feature easy-to-read, die cast beams; heavy-duty welded steel construction; and polystyrene tray.

Model	Capacity	Overall Size	Tray Size
450	130 lb x 1 oz	27"H x 22"W x 21"D	22"L x 14 ³ / ₄ "W x 4 ¹ / ₂ "H
451	65 kg x 20 g	27"H x 22"W x 21"D	22"L x 14 ³ / ₄ "W x 4 ¹ / ₂ "H
459	40 lb x .5 oz/ 17.5 kg x 10 g	27"H x 22"W x 21"D	22"L x 14 ³ / ₄ "W x 4 ¹ / ₂ "H

Model 459-CH same as the 459 but with high-impact, plastic child seat instead of the tray. Seat has comfortable foam pad that is removable for easy cleaning.

Model 243: Capacity: 20 lb x 2 oz
Economical, compact, beam scale
Overall Size: 10¹/₂"H x 21"W x 19"D
Metal Tray: White epoxy - 19"L x 11" W x 2¹/₂"H

Model 2431: Capacity: 10 kg x 50 g

Detecto digital baby scales feature cordless, battery operation, allowing them to be placed where needed. The .7" high contrast, digital display swivels for easy viewing. Additional functions include auto zero, auto shutoff, and weight display lock. Optional 110 VAC adapter available.


Model 8435
Capacity: 130 lb x .1 lb/
40 kg x .02 kg

Tray Size: 22"L x 14³/₄"W x 4¹/₂"H
Overall Size: 20¹/₄"H x 22"W x 22"D

Model 8432-CH
Capacity: 44 lb x .05 lb/
20 kg x .02 kg
Plastic child seat with
removeable pad instead of tray

Model 253
Capacity: 41 lb x 1/4 oz
Sturdy, compact construction; easy-to-read beam
Overall Size: 6³/₄"H x 24"W x 13³/₄"D
Metal Tray: White epoxy - 24"L x 11¹/₂"W x 2³/₄"H

Model 2531
Capacity: 20.5 kg x 10 g



We reserve the right to improve, enhance or modify features and specifications without prior notice.



Sold By:



P.O. Box 151, Webb City, MO 64870 USA
Ph: 417-673-4631 or 1-800-641-2008 Fax: 417-673-5001
Visit our Web site at <http://www.detectoscale.com>

© Copyright 2001 Cardinal Scale Mfg. Co. • Printed in USA • MWC/0701/C162

Bulletin No.C-162